



Packaging Consultants International

1101 State Fair Boulevard, Syracuse, NY 13209

Tel: 315-468-6234

www.pcilabs.com

Fax: 315-468-6235



FDA / IMS #36-B-00155

POUR CONTACT TEST (SM 13.043) REPORT

Testing Lab A2LA Certificate #3497.01

Samples Received From: Frontale Packaging - Syracuse, NY

Date Received: 2/5/2024

Purchase Order: 789456

Report Submitted To: Jane Doe

Time Received: 10:00 AM

Test Report ID: 912324

Date Analyzed: 2/14/2024

Authorization ID: 13

Date Reported: 2/16/2024

PCI Analyst ID(s): 15

PCI Sample ID(s): 1233

Client Sample Identification	No. Items Tested/ Sample	Area Tested (cm ²)	Microbiological Residual Counts ¹			Number Exceeding Standard
			NOTE: Results in bold below exceed standard			
			Bacteria (RBC)		Coliform (RCC)	
December 20, 2023						
3oz. Plastic Lid, Yellow Unprinted WO#201286, Machine: Husky 11	10	33	5 =	0	0	1
			2 =	5 + 45 mold		
			2 =	10		
			1 =	40	Listeria = Not Found	
3oz. Plastic Lid, Purple Unprinted WO#201301, Machine: Husky 12	7	33	3 =	0	0	0
			2 =	1		
			2 =	3	Listeria = Not Found	

Standards:

Residual Bacteria Count (RBC): 1 bacteria colony per cm² of surface area tested (100 cm² tested - RBC limit is 50/item);

Residual Coliform Count (RCC): absent (Ex. 0)

Additional microbiological testing standards are determined by the report recipient organization's Quality & Food Safety Programs

¹ Results stated relate only to the pieces tested, for the sample(s) as received

Melissa R...
Melissa R... Quality Manager
melissa@pcilabs.com

Patrick A. Frontale
Patrick A. Frontale, President
patrick@pcilabs.com

References:

Section 13.043, Surface Contact Methods, from Standard Methods for the Examination of Dairy Products, 17th ed., American Public Health Association, 2004.

Section C. Bacterial Standards and Examination of Single-Service Containers and/or Closures, from Appendix J of the Grade "A" Pasteurized Milk Ordinance, U.S. Department of Health and Human Services, Public Health Service, Food and Drug Administration, current revision.

Reference to Test Method from Standard Methods for the Examination of Dairy Products

Laboratory Internal Tracking Information

Results Section

1 lid exceeded the standard limit

Standards section: Includes notes regarding standards for results, including footnotes of methods outside the scope of accreditation

Reference Section

Sample Report Only