

Test Report Title

	Packaging Consultants International			
	1101 State Fair Boulevard, Syracuse, NY 13209			
	Tel: 315-468-6234	www.pcilabs.com	Fax: 315-468-6235	
FDA / IMS #36-B-00155	ENVIRONMENTAL PETRI DISH REPORT		Testing Lab A2LA Certificate #3497.01	



Petri Dishes Received From:	Frontale Packaging - Syracuse, NY	Quantity Shipped:	16	Purchase Order:	54321
Report Submitted To:	John Dow	Date Shipped:	4/13/2020	Test Report ID:	999820
		Date Received:	4/16/2020	Analyst ID(s):	17 / 15
		Date Reported:	4/30/2020	Sample ID(s):	3998

Sample ID/ Plate #	Microbiological Counts ¹		
	SMA Plate -	PDA Plate -	
	Total Bacteria CFU	Mold CFU	Yeast CFU
1 . Compressed Air - Line 20	1	0	0
2 . Compressed Air - Line 22	8 + 1 mold	1	0
3 . Compressed Air - Line 27	2	0	0
4 . Compressed Air - Oil Room Compressor	2	5	0
5 . Environmental - Warehouse Aisle 2	4 + 2 mold	9	0
6 . Environmental - Warehouse Aisle 14	1 + 5 mold	8	2
7 . Environmental - Packaging	15	20	0
8 . Environmental - Line 37	0	3	0

Standards:
 Environmental microbiological testing standards are determined by the report recipient organization's Quality & Food Safety Programs

¹ The results stated relate only to the items tested, for the sample(s) as received; CFU = Colony Forming Units

Authorizations:

 Melissa Rice, Laboratory Manager melissa@pcilabs.com	 Patrick A. Frontale, President patrick@pcilabs.com
--	--

References:
 Section 6.020, Standard Plate Count Method, from Standard Methods for the Examination of Dairy Products, 17th ed., American Public Health Association, 2004.
 Section 8.112, Acidified Potato Dextrose (Glucose) Agar, from Standard Methods for the Examination of Dairy Products, 17th ed., American Public Health Association, 2004.

Laboratory Internal Tracking Information

Units for test results are Colony Forming Units per plate (CFU/mL)

NOTE: there are no regulatory standards for environmental air plates

Standards section: Includes notes regarding coliform standards

Authorization section: Laboratory Manager contact for technical questions

Reference Section